



White paper | June 2023

Precision medicine for non-small cell lung cancer (NSCLC):

# Emerging trends in molecular analysis



## About the author

Rania Gaspo, B.Pharm. Ph.D.

Global Therapy Area  
Expert, Oncology

Rania Gaspo earned her B.Pharm and Ph.D. at the University of Montreal, Faculty of Pharmacy. She served as a post-doctoral fellow at Montreal Heart Institute before working as a senior scientist at Merck Research Laboratories. Rania then joined Pfizer where she held positions of increasing responsibilities in medical affairs, clinical development, and medical communications in various therapeutic areas, most notably in oncology. She has also supported the launch of many novel medications, led a global team of experienced medical communications managers, and authored more than 25 peer reviewed publications and 60 scientific communications. Rania is also our Cerba Research publications committee chair.

Of all types of cancer, lung cancer caused more deaths worldwide than any other type of cancer.<sup>1</sup> It also causes more deaths than breast and colorectal cancers combined.<sup>1</sup> Patients with the most common lung cancer, non-small cell lung cancer (NSCLC), have higher survival rates than patients with small-cell lung cancer (SCLC), but the outcomes for both remain bleak. According to the National Cancer Institute, the five-year survival rate for NSCLC between 2004 and 2010 was 20.7% compared to only 6.3% for SCLC.<sup>2</sup>

NSCLC patients also have more promising and more precise, treatment options today compared to even 10 years ago. Advances in biomarker and precision medicine have led to the development of immunotherapies and targeted novel treatments that have the potential to improve patient outcomes. In the United States, the FDA has approved dozens of biomarker-driven therapies for NSCLC, including ALK, EGFR, and ROS1 inhibitors.<sup>3</sup> Meanwhile, more than 1,200 treatments are in development.<sup>4</sup>

With precision medicine becoming a more integral component of NSCLC treatment, molecular testing is an important step to help researchers, pathologists, and oncologists understand the genetic underpinnings of this disease. Pathologists commonly use immunohistochemistry (IHC) and/or next-generation sequencing (NGS) to characterize lesions and refine diagnoses, while circulating tumor DNA (ctDNA) analysis has emerged as a promising, non-invasive companion diagnostic.

This white paper explores the roles IHC, NGS, and ctDNA play in NSCLC precision medicine therapy development, as well as how they help healthcare providers determine the most effective treatment regimens for their patients.

## The need for broad molecular profiling

National Comprehensive Cancer Network (NCCN) Clinical Practice Guidelines recommend biomarker testing for patients with newly diagnosed stage IV NSCLC. Oncologists can also consider biomarker testing for patients with squamous cell histology.<sup>5</sup>

NCCN Guidelines suggest broad molecular profiling of multiple biomarkers to identify patients who may benefit from targeted therapies. One of the first steps is testing for biomarker alterations: *EGFR*, *ALK*, *ROS1*, *BRAF*, *KRAS*, *NTRK1/2/3*, *MET ex 14 skipping*, *RET*, *HER2 (ERBB2)* and PD-L1. Most of these genetic mutations or rearrangements have approved targeted or immunotherapies available and/or clinical trial options.<sup>5</sup>

Understanding programmed death-1/programmed-death-ligand 1 (PD-1/PD-L1) levels is especially integral to treatment planning. First-line treatment strategies exist for patients with both negative and positive PD-1/PD-L1 expression, including combinations of surgical intervention, chemotherapy, radiation, and PD-1/PD-L1 inhibitors.<sup>6</sup>

## Immunohistochemistry: A cost-effective screening method for biomarker detection

Pathologists across the globe use IHC to identify and quantify biomarkers for precision medicine treatment and research. Many cancer centers choose conventional, single marker IHC stains for their simplicity, speed of execution, and cost-effectiveness. However, multiplex IHC offers the following additional benefits for NSCLC precision medicine treatment and research.

- **Patient-centric approach** : View and review multiple biomarkers on one slide - no need for patients to undergo additional invasive, often painful rebiopsies.
- **Verify co-expression and spatial organization** : View multiple targets within a preserved tissue architecture.
- **Immune profiling** : Characterize tumors and identify predictive biomarkers for immunotherapy response.
- **Move from preclinical to clinical** : Use multiplex IHC to validate targets, select patients, and characterize efficacy and response. This creates an iterative feedback loop where the ability to predict responses from IHC improves as healthcare providers incorporate the data associated with patient outcomes.

Selected examples of available lung biomarkers at Cerba Research in IHC



**PD-L1**, PD-1, **ALK**, **c-MET**, **ROS1**, **HER2**, **panTRK**, **MEK**, **RET**, IL-1 $\alpha$ , Nkp46, FoxP3, Treg\*, MDSC\*, NK cells\*, CD163, CD56, Ki-67


**Bold** = Targetable Biomarker, \* = Multiplex Assay

## Preferred IHC panels for NSCLC

For efficient NSCLC analysis, Cerba Research recommends two immuno-oncology panels:

- The checkpoint inhibitor (CKI) multiplex panel, which consists of CD3, CD8, PD-1, PD-L1 and a custom marker of your choice.
- The PD-L1 multiplex panel, which consists of CD68, panCK, and PD-L1 markers.


Both panels contain druggable targets and have been appropriately validated on lung specimen. Then, based on the results of these panels, clinical trialists can design a treatment regimen that's optimally targeted to the patient's makeup and cancer type.



PD-L1 multiplex panel. NSCLC

## Available IHC Immuno-oncology multiplex panels at Cerba Research

- T reg (CD3/CD4/CD8/CD25/FoxP3)
- T reg light (CD3/CD8/FoxP3)
- M1/M2 (CD68/CD163/c-maf/pSTAT1)
- CKI (CD3/CD8/PD-1/PD-L1/Custom)
- PD-L1 panel (CD68/panCK/PD-L1)
- Tumor temp (CD3/CD8/Tumor mask) – Chromogenic or Fluorescent
- TRM (CD3/CD8/CD103/CD69/CD49a)





### Next-generation sequencing: Broad-panel technique (Oncopanel)

Using NGS, pathologists and researchers observe the order of nucleotides in targeted DNA or RNA regions or the entire genome of both germline and cancer cells, depending on the scope of the investigation or study, from biomarker research to identification of actionable mutations. When developing targeted NSCLC therapies, NGS allows researchers to detect important biomarkers with high throughput, scalability, and speed. Cerba Research can typically turn around an NGS broad-panel assay for NSCLC, aka oncopanel, in two to three weeks with a relevant and easy-to-read report. This level of efficiency allows patients to start receiving treatment as fast as possible - critically important when you are studying an advanced disease with a low survival rate.

NCCN guidelines recommend NGS broad-panel testing (e.g. oncopanels) as one of the preferred methods to detect genetic alterations and oncogenic driver mutations. Examples of alterations in NSCLC include *ALK*, *ROS1*, *RET* gene fusions, *MET* exon 14 skipping, and *NTRKT1/2/3* gene fusions. *EGFR* is one example of a well-known driver mutation. Other mutations detected by NGS include *BRAF*, *KRAS*, *HER2* and other exploratory biomarkers such as high-level *MET* amplification. Broad oncopanels are one of the most common NGS testing method for detecting these gene alterations and mutations in clinical practice.<sup>5</sup>

With patients’ genetic information, clinical trialists can recommend effective drugs already available or steer patients to appropriate clinical trial options. Until recently, NGS was reserved for reference laboratories and large cancer centers due to cost. Now that the cost to implement NGS has dropped, more hospitals are bringing the technology in-house.<sup>7</sup> Although generally more expensive and less accessible than IHC, both techniques ideally complement themselves to give a full characterization of the tumor.

### Most commonly deployed techniques for NSCLC biomarkers<sup>5</sup>

#### NGS use case

A physician performs a biopsy to obtain a sample of a patient's lung lesion. A molecular pathologist and his team run the sample through a broad-panel NGS validated assay— either DNA, RNA, or both. The result demonstrates the patient suffering from metastatic NSCLC has an *ALK* rearrangement. The physician may opt to prescribe ceritinib, an *ALK* inhibitor approved for use in metastatic NSCLC.<sup>8</sup>

Biomarker	Most commonly deployed techniques	Additional techniques
EGFR	NGS, RT-PCR	Sanger sequencing, single genre
ALK	NGS, IHC, Liquid Biopsy	FISH (reflex), RT-PCR
ROS1	NGS	FISH (reflex), IHC, RT-PCR
BRAF	NGS, RT-PCR, Sanger sequencing	IHC
KRAS	NGS, RT-PCR, Sanger sequencing	
MET	NGS, RNA-based NGS	
RET	NGS, RNA-based NGS	FISH, RT-PCR
NTRK1/2/3	NGS, RNA-based NGS	FISH, IHC, PCR
EGFR T790M	NGS, Liquid Biopsy	
PD-L1	IHC	
HER2	NGS	Sanger sequencing, targeted PCR



## The future is now : Circulating tumor DNA (liquid biopsy)

One of the newer testing methods in oncology, ctDNA, is attracting significant attention due to its non-invasive nature and practicality for early detection, ongoing monitoring, and rapid treatment.<sup>9,10</sup> ctDNA refers to cancer cell DNA that breaks down and releases into the bloodstream. Analyzing genetic alterations and mutations using a ctDNA approach, aka liquid biopsy, helps reduce the need for an invasive tissue biopsy or rebiopsy.

Aside from the ability to obtain genetic information via a blood draw, ctDNA testing offers the following advantages:

- **Serial sampling** : Use in conjunction with or without tissue samples (matched specimens)
- **Ease of use** : Requires less staff to implement
- **Patient comfort** : Patients can provide a sample during an outpatient visit
- **Minimally Invasive** : Reduces the need for biopsy or rebiopsy
- **Tissue is the issue** : 1/5 patients don't have enough tissue
- **Monitor** : Sometimes used to monitor tumor burden longitudinally

Currently, ctDNA diagnostics do not always provide the specificity and/or sensitivity of IHC or NGS on tissue biopsy; however, researchers are currently working to improve the performance of these techniques to meet or exceed IHC and NGS in both dimensions. In addition, ctDNA diagnostics continues to be expensive, reserving the assay to large cancer centers or central laboratories such as Cerba Research. However, the cost effectiveness and specificity/sensitivity is expected to improve over time, much like NGS for tissue biopsy.


## Cerba Research services for NSCLC

With fully equipped histopathology services across the U.S., Europe, and Asia, Cerba Research offers the technology and expertise for expedited NSCLC diagnostics. Serving preclinical to clinical research, Cerba Research has over 250 IHC biomarkers/protocols available and a biobank with a large amount of blocks available and growing. We also develop and validate custom biomarkers according to clinical trial needs. We perform both conventional and multiplex IHC, with the ability to detect up to nine biomarkers in one multiplex IHC panel. Additionally, an international network of board-certified pathologists and consultants analyzes IHC results on demand.

Cerba Research also offers global NGS capabilities with multiple comprehensive oncopanels available for analysis. Our capacity for high throughput, with the ability to sequence more than 1,000 whole human genomes in as little as 10-15 days, provides the information needed to bring targeted medicines to patients faster. We are also one of the few diagnostics solutions providers to offer ctDNA analysis. CerbACT Asia, one of our ImmunoOncology Centers of Excellence, offers multiple ctDNA-based assays using the ACTMonitor® platform. Cerba Research Paris offers the EGFR T790M assay on liquid biopsy and many additional NGS oncopanels, such as the Cerba OncoSign (FFPE), OncoSign ctDNA, and OncoSign 600+, according to your clinical trial needs.

Finally, with Cerba Research, you will have access to:

- A worldwide network of over 800 laboratories operating more than 70 technical platforms
- Over 700 clinical pathology experts in several therapeutic areas (oncology, metabolism, endocrinology, toxicology, genetics, autoimmunity, allergology, pregnancy, infectious disease, neurodegenerative disorders, cardiology, hematology) providing access to patients samples, patients recruitment, and ability to build prospective cohorts
- Broad range of 3000 analysis and more than 1100 instruments available
- A dedicated scientific team to support you from end to end of your project
- An international network of board certified pathologists and consultants involved in IHC scoring, QC and more as aligned with regulatory authorities



## Wide Range of Comprehensive Platforms for your NGS Needs

Cerba Research has experience in a variety of NGS applications, including the following:


- **Solid tumors** : Targeted panels for somatic oncology in solid tumors including, but not limited to, our OncoSign, OncoSign ctDNA and OncoSign 600+ broad panel assays. These panels cover alterations with established, emerging and exploratory value across lung, ovarian, breast, colon, melanoma, bladder, GIST, rare tumors and more
- **Liquid tumors** : Targeted panels for somatic oncology in hematological malignancies. These panels cover mutations with established and emerging value for myeloproliferative neoplasms, chronic myelomonosystic leukemia, myelodysplastic syndroms, acute myeloid leukemia (AML) and more
- **Hereditary cancers** : Oncogenetic panels for hereditary cancers, such as the OncoStar, cover mutations within breast, ovarian, prostate, colon, GI, pancreatic, kidney, neuroendocrine tumors and more
- **Rare disorders** : Germline sequencing such as, but not limited, to BRCA1/BRCA2 with more than 100 off the shelf panels available for numerous rare disorders
- **Constitutional genomics** : Whole exome sequencing (WES) for your constitutional genomic assessments
- **Pathogens identification** : Whole genome sequencing (WGS) for various pathogens identification
- **Viral & Microbiome sequencing** : Various sample type may be analyzed, including but not limited to feces
- **Translocations & Fusions** : With RNA-based NGS sequencing, targeted FISH and IHC screening
- **Cell & Gene therapy** : TCR/BCR immune repertoire sequencing, HLA typing and more for your cell & gene therapy trials

### A full range of library preparation platforms

- Perkin Elmer Sciclone® G3 LiquidHandling Workstations
- Agilent Bravo A & B
- Perkin Elmer Zephyr G3
- Hamilton NGS star
- Hamilton Microlab Nimbus

### And a broad range of high throughput sequencing platforms

- Illumina (MiSeqDx, NextSeq500, NextSeq2000, NextSeq 550Dx)
- PacBio
- Ion Torrent
- Nanopore



## Practice Guidelines Aligned with Cerba Research NGS Offerings

### An Example with Non-small Cell Lung Cancer (NSCLC)

According to the NCCN guidelines<sup>1</sup>, "broad-based genomic testing approaches that efficiently utilize limited biopsy tissue while maximizing diagnostic genomic information are most commonly NGS-based" and "broad genomic profiling may be the most informative approach to examining potential mechanisms of resistance". Check out the table below that outlines relevant lung cancer biomarkers mapped against what Cerba Research may offer.

Lung Cancer Biomarkers	Most commonly deployed <sup>1-4</sup>	Additional Assay(s) <sup>1</sup>	Cerba Research NGS <sup>†</sup>	Cerba Research IHC <sup>†*</sup>	Cerba Research FISH <sup>†</sup>
EGFR	NGS, RT-PCR, Sanger sequencing		X	X	
ALK	NGS, IHC, Liquid Biopsy	FISH (reflex), RT-PCR	X	X	X
ROS1	NGS	FISH (reflex), IHC, RT-PCR	X	X	X
BRAF	NGS, RT-PCR, Sanger sequencing	IHC	X	X	
KRAS	NGS, RT-PCR, Sanger sequencing		X	MEK1	
MET	NGS, RNA-based NGS		X	X	X
RET	NGS, RNA-based NGS	FISH, RT-PCR	X	X	X
NTRK1/2/3	NGS, RNA-based NGS	FISH, IHC, PCR	X	X	X
EGFR T790M	NGS, Liquid Biopsy		X		
PD-L1	IHC			X	
HER2	NGS	Sanger sequencing, targeted PCR	X	X	X

1. NCCN guidelines 2023; 2. Bebb *et al. Curr Oncol* 2021; 3. Cabillic *et al. ESMO Open* 2018;3(6):e419; 4. Li *et al. J Nat Cancer Center* 2021; <sup>†</sup>Cerba Research Data In-house mostly available through the ACTOnco<sup>®</sup>/Cerba Paris (NGS) or CR Montpellier/NY (IHC) or Cerba Paris (FISH); <sup>\*</sup>Validation level may vary; IHC=ImmunoHistoChemistry; NGS=Next-Generation Sequencing

## Discover our Cerba OncoSign FFPE & ctDNA For Solid Tumors\*

Cerba Research OncoSign FFPE & ctDNA panels, alongside our OncoSign600+, covers mutations with established, emerging and exploratory value across lung, ovarian, breast, colon, melanoma, bladder, GIST, rare tumors and more. Our OncoSign FFPE is performed for routine practice in France, in parallel with HRD status which is CE-IVD marked.

### DNA

- AKT1
- ALK
- AR
- BRCA1
- BRCA2
- BAP1
- BRAF
- CDK4
- CDK6
- CTNNB1
- EGFR
- ERBB2
- ESR1
- EIF1AX
- FGFR2
- FGFR3
- FOXL2
- GNA11
- GNAQ
- GNAS
- ALK
- BRAF
- EGFR
- ERBB2
- FGFR1
- FGFR2
- MAP2K1
- MET
- MYD88
- NRAS
- FGFR3
- KRAS
- MET
- NRG1
- NTRK1
- NTRK2
- RAF1
- RET
- SF3B1
- STK11
- NTRK3
- PIK3CA
- PPARG
- ROS1
- RET

### RNA – Fusion

- ALK
- BRAF
- EGFR
- ERBB2
- FGFR1
- FGFR2
- FGFR3
- KRAS
- MET
- NRG1
- NTRK1
- NTRK2
- NTRK3
- PIK3CA
- PPARG
- ROS1
- RET

### Microsatellites

- BAT25
- BAT26
- D2S123
- D5S346
- D17S250
- NR-21
- NR-24
- MONO-27
- ACVR2A
- BTBD7
- DIDO1
- MRE11
- RYR3
- SEC31A
- SULF2

\*The gene list is not exhaustive - please contact us for further informations

## About Sample requirements...

Tumor cell content	>Min 10% and if less then 20% microdissection required
Minimum surface for small samples	4mm <sup>2</sup> with 100% tumor cells
Lymphocyte invasion	<25%
Conditions	Age of block (<1 year)
Fixation time	Ideal: <24h
Type	Ideal: FFPE block or FFPE curls (5 um, n=5) Accepted (not preferred): 10 slides

## Conclusion

As researchers identify more druggable targets through highly specific, sensitive testing methods such as IHC, NGS, and ctDNA-based panels, we expect the development of targeted therapies for lung cancer to progress rapidly and for patient care to improve incrementally. As the process quickens, patients living with this devastating disease will have greater potential to live higher quality and hopefully progression-free lives.

## A Cerba Research snapshot in NSCLC precision medicine



### DNA

- NGS
- Oncopanel
- Custom panels
- ctDNA
- DNA extraction
- Streck tubes
- ddPCR
- qPCR
- Whole exome
- Whole genome



### Tissue

- Multiplex/Simple IHC
- 250+ biomarkers/protocols centralized pathology reading
- Large biobank
- FISH, ISH protocols
- Strong immuno-onco simplex & multiplex panels
- Full histopath service
- Spatial analysis in the tumor microenvironment
- Nanostring™



### RNA

- RNAseq (fusion genes)
- NGS
- Oncopanel
- Custom panels
- rtPCR
- PaxGene®
- Nanostring partnership



### Cell

- FCM panels
- NGF Cytek Aurora under dev (up to 40 colors)
- Immunophenotyping
- Lymphocyte infiltration
- Marker analysis (cell surface/cytoplasmic)




### Protein

- Electrophoresis
- Multiplex cytokine profiling (37-plex)
- 50+ ligand binding assays (ELLA, MSD, ELISPOT, ELISA)

## About Cerba Research

Cerba Research, a strategic provider of diagnostic solutions, supports drug development by leveraging patient data and scientific insight to optimize R&D and commercialization globally. Providing early phase research, clinical development through central laboratory and diagnostic testing, assay and biomarker development and validation – through our global network of specialty laboratories. We partner with government agencies, non-government organizations, as well as pharma and biotech organizations to change the shape of clinical development.

Cerba Research is part of Cerba HealthCare, a leading player in medical diagnosis.





## References

1. World Health Organization. Cancer, Key Facts. Available at <https://www.who.int/news-room/factsheets/detail/cancer>. Accessed: May 7, 2023
2. Trends in Lung Cancer Morbidity and Mortality. American Lung Association, Epidemiology and Statistics Unit, Research and Program Services Division, Available at <https://www.lung.org/lung-health-diseases/lung-disease-lookup/lung-cancer/resource-library/lungcancer-factsheet>. Accessed: May 7, 2023
3. Drugs Approved for Lung Cancer. National Cancer Institute, Available at <https://www.cancer.gov/aboutcancer/treatment/drugs/lung>. Accessed: May 7, 2023
4. GlobalData. Decode the Future of the Healthcare Industry. Available at <https://www.globaldata.com/industries-we-cover/healthcare/>. Accessed: May 7, 2023
5. Non-Small Cell Lung Cancer NCCN Guidelines. National Comprehensive Cancer Network. Available at [https://www.nccn.org/guidelines/category\\_1](https://www.nccn.org/guidelines/category_1). Accessed: May 7, 2023
6. Peng L et al. First-Line Treatment Options for PD-L1-Negative Non-Small Cell Lung Cancer: A Bayesian Network Meta-Analysis. *Front Oncol.* 2021 Jun 23;11:657545. doi: [10.3389/fonc.2021.657545](https://doi.org/10.3389/fonc.2021.657545)
7. Quagliata L. Why Community-Based Hospitals Are Bringing Next-Generation Sequencing In-House. Available at <https://www.clinicallab.com/whycommunity-based-hospitals-are-bringing-nextgeneration-sequencing-in-house-243>. Accessed: May 7, 2023
8. Zykadia® (ceritinib) Highlights of Prescribing Information Revised 3/2019. Available at [https://www.accessdata.fda.gov/drugsatfda\\_docs/label/2019/211225s000lbl.pdf](https://www.accessdata.fda.gov/drugsatfda_docs/label/2019/211225s000lbl.pdf). Accessed: May 7, 2023
9. Fernandez-Lazaro D et al. Liquid Biopsy as Novel Tool in Precision Medicine: Origins, Properties, Identification and Clinical Perspective of Cancer's Biomarkers. *Diagnostics (Basel)*. 2020;10(4):215. doi:[10.3390/diagnostics10040215](https://doi.org/10.3390/diagnostics10040215)
10. FDA Approves First Liquid Biopsy Next-Generation Sequencing Companion Diagnostic Test. FDA News Release, August 7, 2020. Available at <https://www.fda.gov/news-events/press-announcements/fdaapproves-first-liquid-biopsy-next-generationsequencing-companion-diagnostic-test>. Accessed: May 7, 2023

Ghent | Johannesburg | Montpellier | New York | Paris | Shanghai | Ryswick | Schajk | Sydney | Taipei | Tokyo | Rotterdam

